

Wheel chair

Model : 203BJ

	Parameter
Length	101cm
Width	53cm
Height	118~130cm
Package (L*W*H/pcs)	82*61*64
Back Height	54cm
Seat Depth	42cm
Seat Width	46cm
Seat Height from the ground	53cm
Front wheel diameter	20cm
Rear wheel diameter	30cm
Load-bearing	150kg
Net weight	35.5kg



Model : 204BJ

	Parameter	
Length	117cm	
Width	65.5cm	
Height	122~134cm	
Package (L*W*H/pcs)	88*67*84	
Back Height	54cm	
Seat Depth	42cm	
Seat Width	46cm	
Seat Height from the ground	58cm	
Front wheel diameter	20cm	
Rear wheel diameter	62cm	
Load-bearing	150kg	
Net weight	37.7kg	