



OLV-5S / OLV-10S

Home oxygen system

A woman with brown hair, wearing a beige sweater, is smiling and adjusting a white oxygen concentrator. An elderly man with white hair and a beard, wearing a blue button-down shirt and khaki pants, is sitting in a brown armchair and smiling back at her. He is wearing a clear oxygen tube in his nose. The concentrator has a digital display showing "00:00" and a flow control knob. The background shows a cozy living room with framed photos on the wall and a potted plant.

Oxygen
Concentrator

Top 5 Reasons Distributors Trust Us



1,000,000+

Units Sold

Proven success in
100+ countries.

98%

Positive Feedback

High satisfaction &
strong repurchase.

Certified

For Global Markets

FDA, CE, ISO13485, MDSAP

ODM

15+ R&D

Experts driving innovation

30+

Countries

With legal registrations

**Let's grow your
oxygen business
— together.**


Smart Design, Built for Real-Life Use

— Every Detail Focused on Comfort, Safety & Ease —



Large LED display

— simple and clear operation




On/Off Button

— one-touch operation




Humidifier bottle

— gentle, moist oxygen



Medical flowmeter

— accurate flow control



HEPA filter

— blocks bacteria and dust
for safer breathing



Wheels with brakes

— easy to move, very safe

Specification:

Flow Range:	0-5L/min(OLV-5S); 0-10L/min (OLV-10S)
Oxygen Concentration:	90%-96%
Pressure:	20-60kpa
Sound Level:	≤45db (OLV-5S); ≤55db(OLV-10S)
Input Power:	500VA(OLV-5S); 800VA(OLV-10S)
Time Setting:	60-180min
Machine Dimensions:	35*25*67cm (13.7*9.8*26 in)
Warranty:	2 years
Alarm:	Power Failure Alarm Low Voltage Alarm
Package:	1PC/1Carton
G.W/N.W:	25/21kg(OLV-5S); 28/25kg(OLV-10S)
Carton Size:	43*38*75cm
20GP:	234 sets
40GP:	486 sets

Zhengzhou Olive Electronic Technology Co.,Ltd.

Add: 5th Floor, Block D., Building 18,Henan International University ,
Science Park,(East District), 450001 Zhengzhou, Henan, China
Email: Info@zzolive.com
What's app:+86 18638288303
www.zzolive.com

