

PORTABLE OXYGEN CONCENTRATOR



High Oxygen Purity up to: 90%-96%



Light Weight: 1.98kg (4.37lbs with battery)



Pulse Flow 1-6 gears equivalent to 5L/min continuous flow



Battery up to: 1.3-4H±15%



Detachable molecular sieve and battery.



Multiple alarm functions. Lower purity alarm.



**Zhengzhou Olive Electronic
Technology Co.,Ltd.**

www.zzolive.com

E-mail: Info@zzolive.com

What's app: +86 18638288303

Product Model:	OLV-POC5
Charging Voltage:	AC100-240V50/60Hz DC19-21V
Power:	130W
Oxygen Purity:	90%-96% (1st-6th gears)
Flow Rate:	1-6 gears
Product Dimensions:	19.9×18.3×8.6 cm (7.83×7.2×3.39in)
Net Wight:	1.98kg (4.37lbs with battery)
Battery Operation Time:	1.3-4H±15%
Battery Charging Time:	Within 2.5Hours
Operation Sound:	60dB(A)
Standard parts: Pulse machine; Power adapter; Charging cable,Portable bag; 2 sets Oxygen cannula; 2 sets First filtes, 2 sets Hepa filters; User manual.	

