



# E4

## VITAL SIGNS MONITOR



## Performance parameter

- The overall medical white color is novel and generous
- Physical button (machine side switch button, one-button pressure measurement button)
- 4 inch TP fit touch screen
- Sound and light alarm, more convenient for doctors to observe the patient's status
- Built-in 3.7V 2000mAh polymer lithium battery
- Humanized charging base
- Support contact and TYPE-C charging
- Waterproof rating: IPX2
- Host size: 155\*72.5\*28.6

Model \ Parameter	SPO2	ETCO2	NIBP
E4A	✓		
E4B	✓	✓	
E4C	✓		✓
E4D			✓



## Product Features

- 4-inch TFT color high-definition touch screen display.
- Mobile charging , rechargeable treasure, car power charging.
- Support continuous detection of blood oxygen saturation for more than 5 hours.
- Equipped with a support base, which can be charged and stored.
- Provide PI parameters to detect heart disease.
- Real-time data storage, viewing historical data and trend charts.
- Support 500 groups of data storage.
- Small size, easy to carry.
- TP laminated touch screen.
- Dustproof and waterproof grade: IP22.
- Humanized charging base, two-in-one storage and charging.

**Audible alarm**  
It is convenient for medical staff to observe the status of patients in real time

**ETCO2 Sensor Jack  
SpO2 Sensor Jack**  
Multi-parameter monitoring technology, optional CO2, Provide SPO2, Pulse Rate, NIBP monitoring, with excellent parameter monitoring accuracy and anti-interference ability

**Touch screen**  
Standard with touch screen, more convenient for operation  
Support multi-language operation interface

**Charging Basket**

**Charging Base**