



UPGRADE YOUR DEVICE TODAY

COMPARISON OF OUR CURRENT SOLUTIONS VS. LEGACY DEVICES



FEATURE	SPOT VITAL SIGNS 4200 DEVICE <i>(END OF MANUFACTURE)</i>	SPOT VITAL SIGNS LXi DEVICE <i>(END OF SERVICE)</i>	SPOT VITAL SIGNS 4400 DEVICE	CONNEX SPOT MONITOR
Workflows	Spot	Spot	Spot	Spot, Office, and Interval Profiles
BP Averaging	No	No	Yes	Yes
Parameters	Step BP, Temp, SPO ₂	SureBP , Temp, SPO ₂	SureBP , Temp, SPO ₂	SureBP , Temp, SPO ₂ , RRp
SPO ₂ Options	Nellcor, Masimo	Nellcor, Masimo	Nonin	Nellcor, Masimo, Nonin
Temperature Options	SureTemp	SureTemp Plus, Braun PRO 4000	SureTemp Plus	SureTemp Plus, Braun PRO 6000
Custom Parameters/Modifiers	None	Basic (height, weight, pain, RR, BMI)	Basic (height, weight, pain, RR, BMI)	Advanced (20 customizable fields)
Early Warning Scores	No	No	No	Yes
Wireless Connectivity	No	WiFi via Accessory	No	WiFi, Bluetooth
Wired Connectivity	Serial Cable	USB, Serial Cable	USB	USB or Ethernet
Barcode Scanner	No	Yes	No	Yes
Scale Connectivity	No	Yes via USB Dongle	No	Yes via USB Dongle
Display	5" Monochrome LCD	5.8" Monochrome LCD	7" Color Touchscreen	7" Color Touchscreen
Memory	1 Record	50 records	50 records	400 records
Battery Type	Li Ion	Li Ion	Li Ion	Li Ion
Charge Time	12 Hours	10 Hours	4 Hours	4 Hours
Run Time	8 Hours	8 Hours	8 Hours	8 Hours or 17 Hours with APM
Mounting Options	Basic Stand, Wall Mount, IV Pole Mount	Basic Stand, Wall Mount	Desktop, Basic Stand, 777 Wall Mount	Desktop, Basic Stand, APM Stand, 777 or GCX Wall Mount

Contact your Baxter representative for additional information.

[Baxter.com](https://www.baxter.com) | Baxter International Inc.

Baxter, Connex, Spot Vital Signs, SureBP, SureTemp Plus, and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries. Nellcor is a trademark of Covidien LP. Masimo is a trademark of Masimo Corporation. Nonin is a trademark of Nonin Medical, Inc. Braun is a trademark of Braun GmbH.